## **Slides for Preschool Model**

Preschool Children Efron Visual Acuity Test



Test No. 1, SAE-1 Demonstration Right and Left Eye

"Shoe" both eyes.



Test No. 2, SAE-40-2 Right and Left Eye 20/40 Acuity

Different Symbol Each Eye. Also available 20/30 Acuity.



Test No. 3, SAE-40-3 Right and Left Eye 20/40 Acuity


Different Symbol Each Eye. Also available 20/30 Acuity.



Test No. 4, SAE-40-4 Right and Left Eye 20/40 Acuity


Different Symbol Each Eye. Also available 20/30 Acuity.

Primary Grade Children



## Test No. 5, SAS-2 Acuity Right and Left Eye

Conventional Tumbling "E" symbol. Child states which way "E" points in each square. Acuity levels of 20/20, 20/30, 20/40 and 20/50.



Test No. 6, SPS-3 Heterophorias

Muscle imbalance-vertical and lateral is tested at predetermined pass/fail levels.

This test is administered first at far, then at near.



Test No. 7, SPCD-1 - Color Perception

Titmus Pseudo-Isochromatic Plates reliably detect gross color anomalies, but do not classify the type of deficiency.



Test No. 8, SA-BRL-4 - Visual Acuity

Results calibrated at 20/20, 20/30, 20/40, 20/50, 20/70, 20/100 and 20/200.

Tests can be administered binocularly or with either eye occluded.

Improvements have been made to test slides to enhance viewing performance. Patent Pending.

