Slides for Driver Education Model



Test No. 1, SA-BRL-40 Tests Visual Acuity

All lines at 20/40. Selection of test lines eliminates memorization.



Tests ability to read traffic signs in acuity ratings of 20/40 and 20/70.



Test No. 2, SA-BRL-4
Tests Visual Acuity

Results calibrated at 20/20, 20/30, 20/40, 20/50, 20/70, 20/100 and 20/200. Tests can be administered binocularly or with either eye occluded.



Test No. 6, SA-BRL-N2 An alternate to Test No. 2

Sometimes people who cannot recognize letters can recognize numbers.



Test No. 3, SCI-1

Tests ability to recognize traffic colors.



Test No. 7, SLP-1

Measures lateral muscle balance.



Test No. 4, SSR-SD

Tests sign recognition and depth perception.



Test No. 8, SVP-1

Measures vertical muscle balance.

Improvements have been made to test slides to enhance viewing performance. Patent Pending.

